



## COMPUTER TRANSFER STATION

**By simply placing the SilverGuard wand into the Computer Transfer Station (CTS), the data collected in the SilverGuard wand transfers to the connected PC. The CTS is versatile, compact, and is available in both 110V and 220V versions.**

- Reads all SilverGuard Wands with BCS software packages
- Downloads to a PC via USB and Serial Port
- High strength PVC housing
- Powered by AC adapter available in 110 or 220 volts
- Visual and audible confirmation
- Operating temperature: 0°to 80°Celsius
- Memory: 128K SRAM, Processor: DS5001@11MHZ
- Length: 5.69in/14.45cm: Width: 3.58in/9.1cm  
Height: 2.87in/7.3cm
- Reads iButtons directly
- Suitable for all environments, including malls, corrections, border patrol, hospitals, property management, etc.